

BeneCast™ Fibreglass

A fibreglass casting tape that is lightweight yet extremely strong, providing an excellent strength to weight ratio and a smooth finish.

BeneCast™ Fibreglass, is a synthetic fibre based casting tape that provides high quality patient care and comfort. By combining a fibreglass substrate with an advanced water-activated polyurethane resin system. BeneCast™ Fibreglass boasts a low tack application with a smooth finish and excellent lamination providing a strong durable cast. The special air permeability properties also allows moisture to evaporate from the skin helping to reduce irritation.



Features & Benefits

- Lightweight and strong - provides excellent strength to weight ratio
- Air permeability - allows evaporation of moisture from the skin
- Short application time
- Smooth finish
- Low profile

Indications For Use

- Orthopaedic and trauma conditions
- Fracture Management
- Soft tissue injuries
- DFU
- Neurological symptoms



Benecare Medical
 Branson Court, Branson Street,
 Manchester, M40 7FP, United Kingdom
 Email: info@benecaremedical.com
 Tel: +44(0)161 273 6789 Fax: +44 (0) 161 273 1400
www.benecaremedical.com



Product Codes, Sizes & Colours

NHS Code	BeneCare Code	Description	Size (x 3.6m)	Colour
EAF325	FCTB2	Fibreglass Casting Tape	5cm (2")	Blue
EAF326	FCTB3	Fibreglass Casting Tape	7.5cm (3")	Blue
EAF327	FCTB4	Fibreglass Casting Tape	10cm (4")	Blue
	FCTB5	Fibreglass Casting Tape	12.5cm (5")	Blue
EAF329	FCTR2	Fibreglass Casting Tape	5cm (2")	Red
EAF330	FCTR3	Fibreglass Casting Tape	7.5cm (3")	Red
EAF332	FCTR4	Fibreglass Casting Tape	10cm (4")	Red
	FCTR5	Fibreglass Casting Tape	12.5cm (5")	Red
EAF336	FCTG2	Fibreglass Casting Tape	5cm (2")	Green
EAF338	FCTG3	Fibreglass Casting Tape	7.5cm (3")	Green
EAF339	FCTG4	Fibreglass Casting Tape	10cm (4")	Green
EAF340	FCTG5	Fibreglass Casting Tape	12.5cm (5")	Green
EAF341	FCTP2	Fibreglass Casting Tape	5cm (2")	Purple
EAF342	FCTP3	Fibreglass Casting Tape	7.5cm (3")	Purple
EAF343	FCTP4	Fibreglass Casting Tape	10cm (4")	Purple
	FCTP5	Fibreglass Casting Tape	12.5cm (5")	Purple
EAF345	FCTPK2	Fibreglass Casting Tape	5cm (2")	Pink
EAF346	FCTPK3	Fibreglass Casting Tape	7.5cm (3")	Pink
EAF347	FCTPK4	Fibreglass Casting Tape	10cm (4")	Pink
	FCTPK5	Fibreglass Casting Tape	12.5cm (5")	Pink
EAF350	FCTW2	Fibreglass Casting Tape	5cm (2")	White
EAF352	FCTW3	Fibreglass Casting Tape	7.5cm (3")	White
EAF354	FCTW4	Fibreglass Casting Tape	10cm (4")	White
	FCTW5	Fibreglass Casting Tape	12.5cm (5")	White

Call, email or fax through your order with product code, details and quantity and we will dispatch your order immediately.

Buy Online.... www.benecaredirect.com

Benecare Medical
 Branson Court, Branson Street,
 Manchester, M40 7FP, United Kingdom
 Email: info@benecaremedical.com
 Tel: +44(0)161 273 6789 Fax: +44 (0) 161 273 1400
www.benecaremedical.com

