

BeneCast™ Select

The BeneCast™ original polyester casting tape, perfect for conventional casting, as it combines strength and varying levels of rigidity.

BeneCast™ Select, polyester casting tape, is our original quality casting material that is used for conventional primary and secondary casting. BeneCast™ Select combines strength and varying levels of rigidity with soft edges for patient comfort. It also provides a lightweight and breathable cast which is easy to apply and conforms to body contours for easier cast application and increased patient wellbeing.



Features & Benefits

Lightweight and strong - provides excellent strength to weight ratio
Air permeability - allows evaporation of moisture from the skin
Short application time
Smooth finish
Low profile

Indications For Use

Orthopaedic and trauma conditions
Fracture management
Soft tissue injuries
DFU
Neurological Symptoms



BeneCare Medical
 Branson Court, Branson Street,
 Manchester, M40 7FP, United Kingdom
 Email: info@benecaremedical.com
 Tel: +44(0)161 273 6789 Fax: +44 (0) 161 273 1400
www.benecaremedical.com



Product Codes, Sizes & Colours

NHS Code	Benecare Code	Description	Size	Colour	Qty
EAF 192	CTB2	Polyester Casting Tape	5cm	Blue	Box of 10
EAF 193	CTB3	Polyester Casting Tape	7.5cm	Blue	Box of 10
EAF 194	CTB4	Polyester Casting Tape	10cm	Blue	Box of 10
EAF 196	CTR2	Polyester Casting Tape	5cm	Red	Box of 10
EAF 197	CTR3	Polyester Casting Tape	7.5cm	Red	Box of 10
EAF 198	CTR4	Polyester Casting Tape	10cm	Red	Box of 10
EAF 200	CTW2	Polyester Casting Tape	5cm	White	Box of 10
EAF 201	CTW3	Polyester Casting Tape	7.5cm	White	Box of 10
EAF 202	CTW4	Polyester Casting Tape	10cm	White	Box of 10
EAF 203	CTW5	Polyester Casting Tape	12.5cm	White	Box of 10
EAF 204	CTG2	Polyester Casting Tape	5cm	Green	Box of 10
EAF 205	CTG3	Polyester Casting Tape	7.5cm	Green	Box of 10
EAF 206	CTG4	Polyester Casting Tape	10cm	Green	Box of 10
EAF 208	CTP2	Polyester Casting Tape	5cm	Purple	Box of 10
EAF 209	CTP3	Polyester Casting Tape	7.5cm	Purple	Box of 10
EAF 210	CTP4	Polyester Casting Tape	10cm	Purple	Box of 10
EAF 213	CTPK2	Polyester Casting Tape	5cm	Pink	Box of 10
EAF 214	CTPK3	Polyester Casting Tape	7.5cm	Pink	Box of 10
EAF 215	CTPK4	Polyester Casting Tape	10cm	Pink	Box of 10
EAF 217	CTBK2	Polyester Casting Tape	5cm	Black	Box of 10
EAF 218	CTBK3	Polyester Casting Tape	7.5cm	Black	Box of 10
EAF 219	CTBK4	Polyester Casting Tape	10cm	Black	Box of 10
EAF 262	CTO2	Polyester Casting Tape	5cm	Orange	Box of 10
EAF 263	CTO3	Polyester Casting Tape	7.5cm	Orange	Box of 10
EAF 264	CTO4	Polyester Casting Tape	10cm	Orange	Box of 10
EAF2572	CTY2	Polyester Casting Tape	5cm	Yellow	Box of 10
EAF2573	CTY3	Polyester Casting Tape	7.5cm	Yellow	Box of 10
EAF2574	CTY4	Polyester Casting Tape	10cm	Yellow	Box of 10

Call, email or fax through your order with product code, details and quantity and we will dispatch your order immediately.

Buy Online.... www.benecaredirect.com

Benecare Medical
Branson Court, Branson Street,
Manchester, M40 7FP, United Kingdom
Email: info@benecaremedical.com
Tel: +44(0)161 273 6789 Fax: +44 (0) 161 273 1400
www.benecaremedical.com

