

BeneFoot™ 'Original' Post Op Shoe

The BeneFoot Post Op Shoe is at the forefront of patient recovery, helping to increase patient post operative patient safety & stability.

The BeneFoot 'Original' Post Operative Medical Shoe is to be used for diabetic treatment, after post-operative procedures and within wound care, to protect the foot in any situation where swollen or bandaged feet make using normal shoes impractical. Being the only shoe available with a fully rigid 'real' sole and made using only the best materials the BeneFoot 'Original' Post Operative Medical Shoe is more a standard shoe than a medical apparatus with patients able to carry on with normal activities. Every batch of shoes are tested for strength and quality and results indicate that they are twice as durable as industry standards.



Features & Benefits

Dual Density 3/16" Insole – Comfortable, allows for insole modification.
Flared TPR Outsole with rocker bottom – Aids in Ambulation, provides better stability
Padded collar & topline – Reduces heel slippage.
Black upper mesh nylon – Better aesthetics, looks more like a shoe.
Reinforced heel counter – For better medial/lateral stability.
Quality construction – 100% Satisfaction guarantee.
PolyPro shank built into outsole – Allows for better flexion control.
Eight distinct sizes – Better fit & more comfort.
Metal stirrups immobilise the ankle securely and comfortably.

Indications For Use

For recovery after post-operative and tissue procedures.
For use with any size leg cast.
To accommodate large dressings or swellings.
Management of diabetic foot.
Any situation where size of foot prohibits normal footwear.



Benecare Medical
Branson Court, Branson Street,
Manchester, M40 7FP, United Kingdom
Email: info@benecaremedical.com
Tel: +44(0)161 273 6789 Fax: +44 (0) 161 273 1400
www.benecaremedical.com



Product Codes, Sizes & Colours

NHS Code	BeneCare Code	Description	Size	Colour	Quantity
ERF001	POPD	Original Post Op Medical Shoe	Paediatric	Black	12 shoes
ERF002	POWS	Original Post Op Medical Shoe	Women's Small	Black	12 shoes
ERF003	POWM	Original Post Op Medical Shoe	Women's Medium	Black	12 shoes
ERF004	POWL	Original Post Op Medical Shoe	Women's Large	Black	12 shoes
ERF005	POMS	Original Post Op Medical Shoe	Men's Small	Black	12 shoes
ERF006	POMM	Original Post Op Medical Shoe	Men's Medium	Black	12 shoes
ERF007	POML	Original Post Op Medical Shoe	Men's Large	Black	12 shoes
ERF008	POMXL	Original Post Op Medical Shoe	Men's X-Large	Black	12 shoes

Call, email or fax through your order with product code, details and quantity and we will dispatch your order immediately.

Buy Online.... www.benecaredirect.com

Benecare Medical
Branson Court, Branson Street,
Manchester, M40 7FP, United Kingdom
Email: info@benecaremedical.com

Tel: +44(0)161 273 6789 Fax: +44 (0) 161 273 1400
www.benecaremedical.com



BENECARE