



BenePlas™ Aqua & Aqua Colours

This material provides a 100% elastic memory and enhanced mould ability properties, perfect for thumb, wrist, hand and elbow splints.

BenePlas™ Aqua + Aqua Colours have a 100% elastic memory, and goes transparent at moulding temperature enhancing mouldability properties. BenePlas™ Aqua is also manufactured in 5 distinct colours with a coated surface helping to avoid accidental sticking but providing enough tack to hold a temporary bond for the 'wrap n tack' splinting technique. When pressed together the surfaces of the plastic will stick however once cooled, the bond can be popped apart. BenePlas™ Aqua due to the 100% memory can be remodelled numerous times.



Features & Benefits

- BenePlas™ Easy has 100% elastic memory.
- Transparent when reaching moulding temperature of 160oF (70oC)
- Provides assistance with accurate placement of splint to avoid pressure areas.
- Pop apart function when cooled
- Permenant bond can be achieved by pressing firmly and rubbing the two surfaces together.
- Upon cooling, the material will return to its original colour.
- Reheat and remould with ease due to the 100% elastic type memory
- Instant bonding with no finger imprint.
- Extreme elasticity (100%) without tearing.
- Perforations provide varying degrees of ventilation and weight.

Recommended For

- Elbow splints or for any circumferential splints.
- Thumb wrist and elbow splints.



Benecare Medical
Branson Court, Branson Street,
Manchester, M40 7FP, United Kingdom
Email: info@benecaremedical.com
Tel: +44(0)161 273 6789 Fax: +44 (0) 161 273 1400
www.benecaremedical.com



Product Codes, Sizes & Colours

NHS Code	BeneCare Code	Description	Size	Colour	Qty
FNS981	TP104	BenePlas™ Aqua Smooth	460X610X1.6mm	White	1 Sheets
FNS959	TP104H	BenePlas™ Aqua Smooth	460X610X1.6mm	White	2 Sheets
FNS983	TP106	BenePlas™ Aqua Micro Perf	460X610X1.6mm	White	1 Sheets
FNS985	TP111	BenePlas™ Aqua Smooth	460x610x2.4mm	White	1 Sheets
FNS990	TP117	BenePlas™ Aqua Smooth	460x610x3.2mm	White	1 Sheets
FNS991	TP177F	BenePlas™ Aqua Smooth	460x610x3.2mm	White	1 Sheets
FNS788	TP118	BenePlas™ Aqua 1% Perf	460x610x3.2mm	White	1 Sheets
FNS789	TP118F	BenePlas™ Aqua 1% Perf	460x610x3.2mm	White	1 Sheets
FNS999	TP219	BenePlas™ Aqua Colours 19% Perf	460x610x3.2mm	Violet	1 Sheets
FNS1000	TP220	BenePlas™ Aqua Colour Smooth	460x610x2.4mm	Violet	1 Sheets
FNS1002	TP222	BenePlas™ Aqua Colour Smooth	460X610X1.6mm	Violet	1 Sheets
FNS1006	TP228	BenePlas™ Aqua Colour Smooth	460x610x2.4mm	Blue	1 Sheets
FNS1008	TP230	BenePlas™ Aqua Colour Smooth	460X610X1.6mm	Blue	1 Sheets

Call, email or fax through your order with product code, details and quantity and we will dispatch your order immediately.

Buy Online.... www.benecaredirect.com

Benecare Medical
 Branson Court, Branson Street,
 Manchester, M40 7FP, United Kingdom
 Email: info@benecaremedical.com
 Tel: +44(0)161 273 6789 Fax: +44 (0) 161 273 1400
www.benecaremedical.com

