



BenePlas™ Choice

A controlled stretch thermoplastic, suitable for use of experienced clinicians; it provides an excellent conformability making it suitable for a wide range of splints.

BenePlas™ Choice is an excellent choice for the more experienced therapist, due to its intricate conformability. Its controlled stretch provides excellent draping and handling characteristics for a perfect fit to the contours of the limb. It has sufficient rigidity to eliminate the need for reinforcement and has high levels of flexibility. BenePlas Choice is available in both smooth and 1% perforation, and can be used for a variety of splinting techniques.



Features & Benefits

- Controlled stretch that makes larger patterns easy to handle.
- Intricate conformability ensuring precise fit.
- Sufficient rigidity eliminates need for reinforcement.
- Flexibility - stands up to high stress applications (weight bearing splints).
- Colour = Vanilla

Recommended For

- All upper and lower extremity splints
- Static and dynamic splints
- Spasticity, elbow and knee splints.
- Hand, finger and wrist splints.
- Fracture bracing and body jackets.



BeneCare Medical
 Branson Court, Branson Street,
 Manchester, M40 7FP, United Kingdom
 Email: info@benecaremedical.com
 Tel: +44(0)161 273 6789 Fax: +44 (0) 161 273 1400
www.benecaremedical.com



Product Codes, Sizes & Colours

| NHS Code | BeneCare Code | Description | Size | Colour | Qty. |
|----------|---------------|--------------------------|---------------|---------|---------|
| FNS967 | TP701 | BenePlas™ Choice Smooth | 460x610x2.4mm | Vanilla | 1 Sheet |
| FNS968 | TP702 | BenePlas™ Choice 1% Perf | 460x610x2.4mm | Vanilla | 1 Sheet |
| FNS955 | TP703 | BenePlas™ Choice Smooth | 460x610x3.2mm | Vanilla | 1 Sheet |
| FNS956 | TP703F | BenePlas™ Choice Smooth | 460x610x3.2mm | Vanilla | 1 Sheet |
| FNS957 | TP704 | BenePlas™ Choice 1% Perf | 460x610x3.2mm | Vanilla | 1 Sheet |
| FNS1013 | TP704F | BenePlas™ Choice 1% Perf | 460x610x3.2mm | Vanilla | 1 Sheet |

Call, email or fax through your order with product code, details and quantity and we will dispatch your order immediately.

Buy Online.... www.benecaredirect.com

Benecare Medical
Branson Court, Branson Street,
Manchester, M40 7FP, United Kingdom
Email: info@benecaremedical.com

Tel: +44(0)161 273 6789 Fax: +44 (0) 161 273 1400
www.benecaremedical.com

