




BenePlas™ Pro Plus

A thermoplastic material that therapists can work aggressively with due to its special formula that provides enhanced drape-ability.

BenePlas™ Pro-Plus is made from a special formula that provides enhanced drape-ability together with excellent conformity and mould-ability due to its rubber characteristic. It is also extremely rigid when cooled and will not require any reinforcing. BenePlas™ Pro-Plus has a light tack when warm and will remain in place when moulding critical contours.



Features & Benefits

- It will self bond to itself when warm.
- Therapist can aggressively work the material.
- Minimum resistance to stretch.
- High conformability.
- Low control.
- High rigidity.
- Colour = Oyster.

Recommended For:

- Functional hand & wrist splints.
- Upper and lower extremity splints.
- Medium and large splints.
- Base for dynamic outrigger splints.
- Fracture bracing
- Knee splints for immobilisation
- Ankle stirrup splints and foot drop splints.
- Body jackets.
- Burn splints.
- "Serial" static splints and spasticity splints.



Benecare Medical
Branson Court, Branson Street,
Manchester, M40 7FP, United Kingdom
Email: info@benecaremedical.com
Tel: +44(0)161 273 6789 Fax: +44 (0) 161 273 1400
www.benecaremedical.com



Product Codes, Sizes & Colours

NHS Code	BeneCare Code	Description	Size	Qty	Colour
FNS961	TP501	BenePlas Pro Plus™ Smooth	460x610x2.4mm	1 Sheet	Oyster
FNS962	TP502	BenePlas Pro Plus™ 1% Perf	460x610x2.4mm	1 Sheet	Oyster
FNS963	TP505	BenePlas Pro Plus™ Smooth	460x610x3.2mm	1 Sheet	Oyster
FNS963	TP505X	BenePlas Pro Plus™ Smooth	460x610x3.2mm	1 Sheet	Oyster
FNS964	TP506X	BenePlas Pro Plus™ 1% Perf	460x610x3.2mm	1 Sheet	Oyster

Call, email or fax through your order with product code, details and quantity and we will dispatch your order immediately.

Buy Online.... www.benecaredirect.com

Benecare Medical
Branson Court, Branson Street,
Manchester, M40 7FP, United Kingdom
Email: info@benecaremedical.com

Tel: +44(0)161 273 6789 Fax: +44 (0) 161 273 1400
www.benecaremedical.com



BENECARE