

BeneCare™ Finger Mallet Splint

A simple and effective finger splint that provides support to the distal joint of the finger while allowing movement of proximal interphalangeal joint.

BeneCare™ Mallet Finger Splints are designed to support distal joint of finger in extension while permitting unrestricted movement of proximal interphalangeal joint. BeneCare™ Mallet Finger Splints are normally used without padding and is taped to middle phalanx, elevating distal joint to properly extended position.



Features & Benefits

- Made from high quality rigid plastic.
- Ventilated design for greater comfort.
- Removeable stay can be shaped.
- 7 sizes available.

Indications For Use

- To support distal joint of finger in extension.
- Permits unrestricted movement of proximal interphalgeal joint.



BeneCare Medical
Branson Court, Branson Street,
Manchester, M40 7FP, United Kingdom
Email: info@benecaremedical.com
Tel: +44(0)161 273 6789 Fax: +44 (0) 161 273 1400
www.benecaremedical.com



Product Codes, Sizes & Colours

NHS Code	BeneCare Code	Description	Size	
GMB243	FMS1	Mallet Finger Splint	Size 1	X-Small
GMB244	FMS2	Mallet Finger Splint	Size 2	Small
GMB245	FMS3	Mallet Finger Splint	Size 3	Medium
GMB246	FMS4	Mallet Finger Splint	Size 4	Medium/Large
GMB247	FMS5	Mallet Finger Splint	Size 5	Large
GMB248	FMS5.5	Mallet Finger Splint	Size 5.5	Large/X-Large
GMB249	FMS6	Mallet Finger Splint	Size 6	X-Large
GMB250	FMS7	Mallet Finger Splint	Size 7	XX-Large

Call, email or fax through your order with product code, details and quantity and we will dispatch your order immediately.

Buy Online.... www.benecaredirect.com

Benecare Medical
 Branson Court, Branson Street,
 Manchester, M40 7FP, United Kingdom
 Email: info@benecaremedical.com
 Tel: +44(0)161 273 6789 Fax: +44 (0) 161 273 1400
www.benecaremedical.com

