


# BenePlas™ Neo Prene Sheets

A lighter alternative to traditional thermoplastic splinting and perfect for creating flexible splints tailored to individual patient contours.

BenePlas™ Neoprene Sheets allow you to create your own flexible splints tailored to your patients needs. BenePlas™ Neoprene has a degree of elasticity enabling controlled movement and consists of a nylon lining on one side and a terry cotton material on the other over a rubber core ensuring high levels of warmth and comfort. BenePlas™ Neoprene can be used as strapping on a splint, or can be used to make soft splints by fusing pieces together with seam tape and iron.



## Features & Benefits

Size of sheet 460 x 610 x 3.2mm
Terry cotton cover
Lighter in weight than thermoplastic splints.

## Indications For Use

Can be used to create flexible splints tailored to patients needs.
Can be used for strapping on a splint.
Can be used for soft splints by fusing pieces together.



Benecare Medical  
Branson Court, Branson Street,  
Manchester, M40 7FP, United Kingdom  
Email: [info@benecaremedical.com](mailto:info@benecaremedical.com)  
Tel: +44(0)161 273 6789 Fax: +44 (0) 161 273 1400  
[www.benecaremedical.com](http://www.benecaremedical.com)



## Product Codes, Sizes & Colours

NHS Code	BeneCare Code	Description	Size	Colour	Qty
	NSBK	BenePlas Neoprene Sheets	41 x 61cm x 3.2mm	Black	1 Sheet
	NSBE	BenePlas Neoprene Sheets	41 x 61cm x 3.2mm	Beige	1 Sheet
	NSBL	BenePlas Neoprene Sheets	41 x 61cm x 3.2mm	Blue	1 Sheet
	NSR	BenePlas Neoprene Sheets	41 x 61cm x 3.2mm	Red	1 Sheet
	NSP	BenePlas Neoprene Sheets	41 x 61cm x 3.2mm	Purple	1 Sheet

Call, email or fax through your order with product code, details and quantity and we will dispatch your order immediately.

**Buy Online.... [www.benecaredirect.com](http://www.benecaredirect.com)**

Benecare Medical  
Branson Court, Branson Street,  
Manchester, M40 7FP, United Kingdom  
Email: [info@benecaremedical.com](mailto:info@benecaremedical.com)

Tel: +44(0)161 273 6789 Fax: +44 (0) 161 273 1400  
[www.benecaremedical.com](http://www.benecaremedical.com)

