

# BeneCare™ Neo-Wrap Brace

The ideal brace for patients suffering with rheumatoid arthritis and people requiring extra support surrounding the knee as it increases medial and lateral support.

The BeneCare™ Neo-Wrap Brace front opening design and easy opening tabs on the Velcro™ straps makes this brace suitable for the arthritic or elderly patient, who may suffer from limited dexterity. The brace incorporates hyper stop hinges provide medial and lateral support for unstable knee, while allowing flexion. The interal access to the hinge allows for easy removal of the hinge for laundering.



## Features & Benefits

Front opening
Universal left and right.
Contrasting proximal & distal strap to counter rotation.
Open popliteal area for added comfort during flexion.
Removable hyper stop hinges.
Easy grip straps for patients with poor dexterity.
Fully adjustable front & rear straps

## Indications For Use

Increasing medial/lateral support
Rheumatoid arthritis
Mild to moderate hyperextension
Mild to moderate valgus/varus deformity.



BeneCare Medical  
 Branson Court, Branson Street,  
 Manchester, M40 7FP, United Kingdom  
 Email: [info@benecaremedical.com](mailto:info@benecaremedical.com)  
 Tel: +44(0)161 273 6789 Fax: +44 (0) 161 273 1400  
[www.benecaremedical.com](http://www.benecaremedical.com)



## Product Codes, Sizes & Colours

NHS Code	BeneCare Code	Description	Size	Colour
GRS2100	SKKBXS	Knee Support - Short	X-Small	Beige
GRS2084	SKKBS	Knee Support - Short	Small	Beige
GRS2085	SKKBM	Knee Support - Short	Medium	Beige
GRS2086	SKKBL	Knee Support - Short	Large	Beige
GRS2098	SKKBXL	Knee Support - Short	X-Large	Beige
GRS2099	SKKBXXL	Knee Support - Short	XX-Large	Beige
GRA424	BSKKBXS	Knee Support - Short	X-Small	Black
GRA521	BSKKBS	Knee Support - Short	Small	Black
GRA520	BSKKBM	Knee Support - Short	Medium	Black
GRA519	BSKKBBL	Knee Support - Short	Large	Black
GRA522	BSKKBXL	Knee Support - Short	X-Large	Black
GRA425	BSKKBXXL	Knee Support - Short	XX-Large	Black

NHS Code	BeneCare Code	Description	Size	Colour
GRA512	LSKKBXS	Knee Support - Long	X-Small	Beige
GRA210	LSKKBS	Knee Support - Long	Small	Beige
GRA509	LSKKBM	Knee Support - Long	Medium	Beige
GRA508	LSKKBBL	Knee Support - Long	Large	Beige
GRA511	LSKKBXL	Knee Support - Long	X-Large	Beige
GRA513	LSKKBXXL	Knee Support - Long	XX-Large	Beige
GRA506	BLSKKBXS	Knee Support - Long	X-Small	Black
GRA500	BLSKKBS	Knee Support - Long	Small	Black
GRA501	BLSKKBM	Knee Support - Long	Medium	Black
GRA502	BLSKKBBL	Knee Support - Long	Large	Black
GRA505	BLSKKBXL	Knee Support - Long	X-Large	Black
GRA507	BLSKKBXXL	Knee Support - Long	XX-Large	Black

Call, email or fax through your order with product code, details and quantity and we will dispatch your order immediately.

**Buy Online.... [www.benecaredirect.com](http://www.benecaredirect.com)**

Benecare Medical  
Branson Court, Branson Street,  
Manchester, M40 7FP, United Kingdom  
Email: [info@benecaremedical.com](mailto:info@benecaremedical.com)  
Tel: +44(0)161 273 6789 Fax: +44 (0) 161 273 1400  
[www.benecaremedical.com](http://www.benecaremedical.com)

