

BeneCare™ Non-Woven Swabs

Super soft swabs with a high absorbancy threshold make them perfect to use for covering and cleaning minor wounds.

BeneCare™ Medical's Non Woven Swabs are made from top grade absorbent fabrics. They have a standard 70% viscose - 30% Polyester fibre mix, with a 30gsm weight. These swabs provide improved strength, absorbency, and scrubbing texture compared to traditional cotton gauze swabs.



Features & Benefits

- Higher absorbency than 12-ply cotton gauze swabs.
- Less linting than traditional gauze.
- Holds its shape and texture better than cotton gauze swabs.
- Reduces snagging on structures and provides excellent wound release.
- Versatile use, cleans, absorbs and wicks.
- Conforms to Drug Tariff Specification 28.
- Soft and gentle yet durable when wet.

Indications For Use

- Clean and cover minor wounds.
- To absorb secretion.
- To treat secondary healing wounds.



Benecare Medical
Branson Court, Branson Street,
Manchester, M40 7FP, United Kingdom
Email: info@benecaremedical.com
Tel: +44(0)161 273 6789 Fax: +44 (0) 161 273 1400
www.benecaremedical.com



Product Codes, Sizes & Colours

NHS Code	BeneCare Code	Size	Qty
	NWS5	5x5cm	500
	NWS7.5	7.5x7.5cm	500
	NWS10	10x10cm	500

Call, email or fax through your order with product code, details and quantity and we will dispatch your order immediately.

Buy Online.... www.benecaredirect.com

Benecare Medical
Branson Court, Branson Street,
Manchester, M40 7FP, United Kingdom
Email: info@benecaremedical.com
Tel: +44(0)161 273 6789 Fax: +44 (0) 161 273 1400
www.benecaremedical.com



BENECARE