


BeneCast™ Pre-Cut Moulded Splints

The pre-cut moulded splints are a fast and efficient product to use when splinting, perfect to use in A&E, theatres and plaster rooms.

BeneCast™ Pre-Cut Moulded Splints are an excellent way of making a rigid splint quickly and effectively. Simply submerge the splint in water for 10 seconds and apply to the body wrapping with a bandage to form its shape. Allow 20-30 minutes for the splint to set and the patient will now be able to move and remove the splint as required.



Features & Benefits

- Can be used quickly and efficiently in A&E, Theatres, Plaster Rooms, Hand Therapy.
- Can be used for ambulance services in emergency situations.
- 5-times lighter and stronger than Plaster Casting.
- Fast and efficient method to make a removable cast.
- Can be secured by either bandage or hook & loop tape.
- Pre-cut to its individual shape that splint requires minimal trimming or preparation.
- Excellent conformability to mould to contours of the limb.
- 100% penetration of X-Ray.

Indications For Use

- Controlled elbow flexion/extension following surgery or trauma.
- Treatment of elbow dislocations.



Benecare Medical
Branson Court, Branson Street,
Manchester, M40 7FP, United Kingdom
Email: info@benecaremedical.com
Tel: +44(0)161 273 6789 Fax: +44 (0) 161 273 1400
www.benecaremedical.com



Product Codes, Sizes & Colours

NHS Code	BeneCare Code	Description	Size	Qty
EAF 300	PCMS001	BeneCast Pre-Cut Splints	Colles	10 Splints
EAF 301	PCMS002	BeneCast Pre-Cut Splints	Volar Splint	10 Splints
EAF 302	PCMS003	BeneCast Pre-Cut Splints	Volar Splint	5 Splints
EAF 303	PCMS004	BeneCast Pre-Cut Splints	Spica Splint	10 Splints
EAF 304	PCMS005	BeneCast Pre-Cut Splints	Long Arm Splint	10 Splints
EAF 305	PCMS006	BeneCast Pre-Cut Splints	Long Leg	5 Splints
EAF 306	PCMS007	BeneCast Pre-Cut Splints	Short Leg Adult	5 Splints
EAF 307	PCMS010	BeneCast Pre-Cut Splints	Short Leg Child	5 Splints

Call, email or fax through your order with product code, details and quantity and we will dispatch your order immediately.

Buy Online.... www.benecaredirect.com

Benecare Medical
Branson Court, Branson Street,
Manchester, M40 7FP, United Kingdom
Email: info@benecaremedical.com
Tel: +44(0)161 273 6789 Fax: +44 (0) 161 273 1400
www.benecaremedical.com



BENECARE