

BeneCare™ Pro Universal Wrist Brace

A lightweight and durable wrist brace perfect to use for tendonitis, healing fractures, soft tissue injuries and sports injuries.

The BeneCare Pro Universal Wrist Brace is constructed of a lightweight and durable fabric facilitating comfort and movement. Added to this the soft Lycra® lined interior is extremely breathable and keeps fresh. The Universal Wrist Brace's ingenious construction allows it to be used injuries on the both the right and the left hand. just slide the shaped bar from one side of the brace and insert it into the other to change hand usage. The Universal Wrist Brace has an inbuilt sock style control system allowing for any user to apply single handedly, helping patients continue with normal life.



Features & Benefits

- Universal Left & Right.
- Extremely breathable Lycra®- lined interior.
- Lightweight durable fabric construction.
- Contoured along distal palmer crease to allow for 90°.
- MCP flexion Precisely sculpted narrow web space prevents pinching.
- Removable / adjustable stay.
- Easy to use contact closure strap.

Indications For Use

- Tendonitis healing fracture.
- Soft tissue injuries



Benecare Medical
Branson Court, Branson Street,
Manchester, M40 7FP, United Kingdom
Email: info@benecaremedical.com
Tel: +44(0)161 273 6789 Fax: +44 (0) 161 273 1400
www.benecaremedical.com



Product Codes, Sizes & Colours

NHS Code	BeneCare Code	Description	Size	Colour
GQV164	UWSS	Universal Wrist Brace	Small	Black
GQV163	UWSM	Universal Wrist Brace	Medium	Black
GQV162	UWSL	Universal Wrist Brace	Large	Black
GQV165	UWSXL	Universal Wrist Brace	X-Large	Black

Call, email or fax through your order with product code, details and quantity and we will dispatch your order immediately.

Buy Online.... www.benecaredirect.com

Benecare Medical
Branson Court, Branson Street,
Manchester, M40 7FP, United Kingdom
Email: info@benecaremedical.com
Tel: +44(0)161 273 6789 Fax: +44 (0) 161 273 1400
www.benecaremedical.com

