

BeneCare™ Silicon Heel Cups

Engineered to help relieve pain from plantar fasciitis, heel spurs, bursitis, Achilles tendon and pressure ulcers.

The BeneCare silicone/gel heel cups use dual durometer for maximum comfort while standing, walking or running, the medium durometer silicone in the main body of the heel cup efficiently absorbs impact and recovers quickly. The silicone heel cups are suitable for indications such as: plantar fasciitis, heel spurs, heel pain, bursitis, Achilles tendonitis and pressure ulcers.



Features & Benefits

Uses dual durometer for maximum comfort while standing, walking or running.
Medium-durometer silicone in main body efficiently absorbs impact and recovers quickly.
Softer silicone blue zone under heel centre
Maximises shock absorption and comfort at point of impact
Retains shape - will not deform or degrade with use
Washable and reusable.

Indications For Use

Plantar fasciitis
Heel Spurs
Heel Pain
Shock absorbent
Bursitis
Achilles tendonitis
Pressure ulcers



Benecare Medical
Branson Court, Branson Street,
Manchester, M40 7FP, United Kingdom
Email: info@benecaremedical.com
Tel: +44(0)161 273 6789 Fax: +44 (0) 161 273 1400
www.benecaremedical.com



Product Codes, Sizes & Colours

NHS Code	BeneCare Code	Size	UK Size	EU Size
ERU4558	SHC-S	Small	3-5	35-38
ERU4559	SHC-M	Medium	6-8	38.5-42.5
ERU4560	SHC-L	Large	9-11	43-45.5
ERU4651	SHC-XL	X-Large	12-14	46-49

Call, email or fax through your order with product code, details and quantity and we will dispatch your order immediately.

Buy Online.... www.benecaredirect.com

Benecare Medical
Branson Court, Branson Street,
Manchester, M40 7FP, United Kingdom
Email: info@benecaremedical.com

Tel: +44(0)161 273 6789 Fax: +44 (0) 161 273 1400
www.benecaremedical.com



BENECARE

