

BeneCare™ Swedish AFO

A low profile leaf spring design to substantially correct flaccid drop foot, through the use of heat and modification conforming to individual patient contours.

The BeneCare Swedish AFO is a comfortable, low profile prefabricated leaf spring design to correct flaccid foot drop, especially on post-stroke patients. Injection moulded polypropylene allows for different thicknesses throughout the orthosis, thicker on the vertical aspect for rigidity, thinner footplate for ease of trimming, comfort and flexibility.



Features & Benefits

- Injection-moulded for a high quality finish.
- Easily trimmed footplate may be heat modified for minor fitting adjustments.

Indications For Use

- Drop foot.
- Mild drop foot secondary to other neurological pathologies



Benecare Medical
Branson Court, Branson Street,
Manchester, M40 7FP, United Kingdom
Email: info@benecaremedical.com
Tel: +44(0)161 273 6789 Fax: +44 (0) 161 273 1400
www.benecaremedical.com



Product Codes, Sizes & Colours

NHS Code	BeneCare Code	Size	UK	EU
GRD1327	SAFORS	Small Right	1-3.5	33-36
GRD1328	SAFORM	Medium Right	4-7	37-40
GRD1329	SAFORL	Large Right	7.5-9.5	41-44
GRD1330	SAFORXL	X-Large Right	10-14	45-48
GRD1323	SAFOLS	Small Left	1-3.5	33-36
GRD1324	SAFOLM	Medium Left	4-7	37-40
GRD1325	SAFOLL	Large Left	7.5-9.5	41-44
GRD1326	SAFOLXL	X-Large Left	10-14	45-48

Call, email or fax through your order with product code, details and quantity and we will dispatch your order immediately.

Buy Online.... www.benecaredirect.com

Benecare Medical
Branson Court, Branson Street,
Manchester, M40 7FP, United Kingdom
Email: info@benecaremedical.com
Tel: +44(0)161 273 6789 Fax: +44 (0) 161 273 1400
www.benecaremedical.com



BENECARE