

BeneCare™ ECO Walker Boot

An ultra light weight walker boot that is impact resistant and offers a secure fit while enhancing protection and comfort to patient.

The BeneFoot™ ECO Walker Boot is an ultra light-weight walker boot providing support and comfort to patients requiring short term immobilisation. The ECO Walker Boot comprises a rocker sole removable foam lining helping improve short term, rehabilitation capabilities.



Features & Benefits

Hook & Loop closure provides easy adjustment
Ultra low profile design eliminates height difference and facilitates natural gait
Fixed stirrups immobilise the ankle securely and comfortably
Rocker bottom sole allows easier ambulation
Padded soft foam keeps lower leg comfortable
Ultra lightweight design makes it suitable for short term immobilisation

Recommended For

Replaces the need for a heavy cast
Stable fracture of fibula, calcaneus and malleolus
Post ligament, soft tissue and tendon surgery
Post-surgery management

NHS Code	BeneCare Code	Description	Size	Colour	Qty
ERE256	FWS	Fixed Walker Eco	Small	Black	10 Boots
ERE257	FWM	Fixed Walker Eco	Medium	Black	10 Boots
ERE258	FWL	Fixed Walker Eco	Large	Black	10 Boots



Benecare Medical
Branson Court, Branson Street,
Manchester, M40 7FP, United Kingdom
Email: info@benecaremedical.com
Tel: +44(0)161 273 6789 Fax: +44 (0) 161 273 1400
www.benecaremedical.com

