

BeneCast™ Flex

A semi rigid polyester based casting tape, that provides semi-rigid support to the limb, while allowing some soft tissue function.

BeneCast™ Flex, semi-rigid casting tape, is made out of polyester based materials which make BeneCast™ Flex stronger and more durable than other semi-rigid casting tapes. The BeneCast™ Flex never fully hardens yet remains resilient, even in its cured state. It provides a semi-rigid support of the limb, while allowing some degree of soft tissue function, and will revert to its cured shape after deformation. Made from multi-directional stretch bandage that conforms to the body contours without folds or tucks. It enables easy application with soft, smooth edges for patient comfort. BeneCast™ Flex Semi Rigid is easy to unroll, can be applied either wet or dry, and has tack free additives to increase control and reduce application time.



Features & Benefits

- BeneCast™ Flex is easy to use.
- Easy to unroll and apply without the need for folds and tucks.
- BeneCast™ Flex is tack-free and clean to apply.
- High porosity and good moisture resistance allows Flex™ casts to breathe and dry out easily.
- Reduces risk of body odour and skin maceration under the cast.
- Easy to remove either by unwrapping or cutting with bandage scissors .
- Apply wet or dry.
- Ideal for paediatric patients.
- Flexible stabilisation.
- Low profile & quick setting.

Indications For Use

- Orthopaedic and trauma conditions
- Fracture management.
- Soft tissue injuries
- DFU
- Neurological symptoms



BeneCare Medical
 Branson Court, Branson Street,
 Manchester, M40 7FP, United Kingdom
 Email: info@benecaremedical.com
 Tel: +44(0)161 273 6789 Fax: +44 (0) 161 273 1400
www.benecaremedical.com



Product Codes, Sizes & Colours

NHS Code	BeneCare Code	Description	Size	Colour	Qty
EAF 270	CTSRB2	Flexible Casting Tape	5cm (2")	Blue	Box of 10
EAF 274	CTSRB3	Flexible Casting Tape	7.5cm (3")	Blue	Box of 10
EAF 278	CTSRB4	Flexible Casting Tape	10cm (4")	Blue	Box of 10
	CTSRR2	Flexible Casting Tape	5cm (2")	Red	Box of 10
EAF 2578	CTSRR3	Flexible Casting Tape	7.5cm (3")	Red	Box of 10
	CTSRR4	Flexible Casting Tape	10cm (4")	Red	Box of 10
EAF 2579	CTSRY3	Flexible Casting Tape	7.5cm (3")	Yellow	Box of 10
	CTSRBK2	Flexible Casting Tape	5cm (2")	Black	Box of 10
EAF 2576	CTSRBK3	Flexible Casting Tape	7.5cm (3")	Black	Box of 10
EAF 2577	CTSRG3	Flexible Casting Tape	7.5cm (3")	Green	Box of 10
EAF 271	CTSRP2	Flexible Casting Tape	5cm (2")	Purple	Box of 10
EAF 275	CTSRP3	Flexible Casting Tape	7.5cm (3")	Purple	Box of 10
EAF 279	CTSRP4	Flexible Casting Tape	10cm (4")	Purple	Box of 10
EAF 272	CTSRPK2	Flexible Casting Tape	5cm (2")	Pink	Box of 10
EAF 276	CTSRPK3	Flexible Casting Tape	7.5cm (3")	Pink	Box of 10
EAF 280	CTSRPK4	Flexible Casting Tape	10cm (4")	Pink	Box of 10
EAF 273	CTSRW2	Flexible Casting Tape	5cm (2")	White	Box of 10
EAF 277	CTSRW3	Flexible Casting Tape	7.5cm (3")	White	Box of 10
EAF 281	CTSRW4	Flexible Casting Tape	10cm (4")	White	Box of 10

Call, email or fax through your order with product code, details and quantity and we will dispatch your order immediately.

Buy Online.... www.benecaredirect.com

Benecare Medical
 Branson Court, Branson Street,
 Manchester, M40 7FP, United Kingdom
 Email: info@benecaremedical.com
 Tel: +44(0)161 273 6789 Fax: +44 (0) 161 273 1400
www.benecaremedical.com

