

BeneCast™ Plaster Of Paris

BeneCast™ POP provides excellent conformability and accommodates the needs of busy fracture clinics, A&E and theatre plaster rooms.

BeneCast™ POP combines all the technician/patient benefits of the finest leno-weave gauze, impregnated with a specially formulated blend of alpha and beta Plaster of Paris crystals incorporating high-grade binders and accelerants. It is then spooled onto a rigid plastic core of layered into slabs.



Features & Benefits

- Excellent wetting-out of the bandage upon immersion into water.
- Full set cast after 3-5 minutes (usual weight bearing advice applies).
- Excellent moulding to the body contours.
- Strong and durable.
- Leno-weave has excellent conformance. Reducing distortion and creasing during application.
- Rigid structure BeneCast™ sets to a hard, rigid, aesthetically please structure.
- Perfect for busy Fracture Clinics, A&E and Theatre Plaster Rooms where rapid immobilisation of injuries is vital.
- Leno-weave fabric ensures extra stability during wet stage, making it ideal for moulding around awkward contours.

Indications For Use

- Orthopaedic and trauma conditions
- Fracture management
- Soft tissue injuries
- DFU



Benecare Medical
 Branson Court, Branson Street,
 Manchester, M40 7FP, United Kingdom
 Email: info@benecaremedical.com
 Tel: +44(0)161 273 6789 Fax: +44 (0) 161 273 1400
www.benecaremedical.com



Product Codes, Sizes & Colours

NHS Code	BeneCare Code	Description	Size	Qty
EAA 020	P.O.P.2	Plaster Of Paris - Rolls	5cm x 2.7m	Box of 24
EAA 022	P.O.P.3	Plaster Of Paris - Rolls	7.5cm x 2.7m	Box of 24
EAA 035	P.O.P.4	Plaster Of Paris - Rolls	10cm x 2.7m	Box of 12
EAA 036	P.O.P.6	Plaster Of Paris - Rolls	15cm x 2.7m	Box of 12
EAA 037	P.O.P.8	Plaster Of Paris - Rolls	20cm x 2.7m	Box of 12
EAA 038	P.O.P.S4	Plaster Of Paris - Slabs	10cm x 20m	Single
EAA 039	P.O.P.S6	Plaster Of Paris - Slabs	15cm x 20m	Single
EAA 048	P.O.P.S8	Plaster Of Paris - Slabs	20cm x 20m	Single

Call, email or fax through your order with product code, details and quantity and we will dispatch your order immediately.

Buy Online.... www.benecaredirect.com

Benecare Medical
 Branson Court, Branson Street,
 Manchester, M40 7FP, United Kingdom
 Email: info@benecaremedical.com
 Tel: +44(0)161 273 6789 Fax: +44 (0) 161 273 1400
www.benecaremedical.com

